

News from Jinan TE _ 201709


Ipad Controller for Testing Machine

After abundant trail running, **TE MaxTest software is proven to be fully compatible with Ipad controller for testing machines of Jinan TE brand**, including the electromechanical motor driven UTM WDW series, servo hydraulic UTM WAW series & HLC series, the TYE series & CTEN series compression tester, and the HTT horizontal tensile testing system etc.


Some advantage over the traditional digital controller is listed below;

1. All the tests are display in real time, including the test data and the curves;
2. The controller with 32 bits ARM provides higher arithmetic speed and data processing speed;
3. WIFI is available, and it also supports USB2.0, Ethernet, and RS485;
4. Data can be stored as TXT for analysis;
5. Big memory space for more test result.
6. More elegant appearance, easy & simple operation;

This new Ipad controller system will be supplied with our machines since this October 2017, and the previous digital controller will be stopping production.



Ipad with the software interface



Software display



Testing machines compatible with the Ipad controller