

QC-528K Universal Testing Machines (500N)

Standard

Description

This model is designed to be available to mount with different brand/type of force gauge, test stand is provided with displacement display, test speed with digital setting, it is the best choice for small specimen test. Light weight stand, small size without taking space.

Industry

Plastic, Paper, Textile, Tape, Packaging

Appearance



Specification

Model	QC-528K
Max. capacity	500 N
Test Stroke	400 mm(w/o grip)
Stroke resolution	0.05 mm
Speed range	6 ~ 600 mm/min (Digital setting)
Hardware safety protection	Up and down limit switch, emergency stop button protection
Machine power	Servo motor
Power	Single phase 110~240VAC · 5A
Dimension	45 × 27 × 88 cm
Weight	28 kg

Related Products

1. QC-508A2 Computerized material testing machine
2. QC-513A2 Computerized material testing machine
3. QC-508M1 Computerized material testing machine
4. QC-513M1 Computerized material testing machine
5. QC-508E Computerized material testing machine
6. All kind of grips